

# WaveX.Net Software Library



## 1. Introduction

The WaveX.NET software library implements all the functionalities required to fully control WaveX and Waveplus devices through a Personal Computer. The communication with the WaveX and Waveplus devices happens by USB cable, following the USB 3.0 or USB 2.0 specifications, by using a kernel mode USB function driver, which is Windows Driver Foundation (WDF) compliant. The WaveX.NET library can be used in C# or C++ Visual Studio projects for .NET platforms with Windows 11/10 (32/64 bit).

## 2. General Description

WaveX.NET is a multi-threading library developed in C# language which allows:

- The WaveX and Waveplus devices configuration:
  - data acquisition process configuration
  - sensors configuration
  - analog outputs configuration (available for WaveX system)
- Continuous data acquisition and acquired data integrity evaluation
- Monitoring of external Start/Stop triggers
- Impedance check on the EMG electrodes (available for Waveplus system)
- Offset compensation of signals produced by the accelerometers
- Inertial sensors calibration


Being thread-safe, WaveX.NET realizes also a state machine in order to guarantee higher software reliability, avoiding wrong sequences in the access to the library functionalities, since it generates [exceptions](#) in case of errors occurred during the implementation of its own functionalities.

The WaveX.NET includes the following assemblies:

*WaveX.dll*  
*WaveX.Common.dll*  
*WaveX.BSys.dll*  
*WaveX.PSys.dll*  
*WaveX.XSys.dll*  
*PicoBlue.DaqSys.dll*  
*CyUSB.dll*

## 3. Architecture

Scheme 1 shows the main software modules.



### 3.1 The DaqSystem Module

The DaqSystem module (namespace [WaveX](#)) includes the homonymous class whose methods and properties allow full control of the WaveX and Waveplus devices. This class implements the IDaqSystem interface (namespace [WaveX.Interfaces](#)) which includes the following properties, methods, and event handlers:

Red Rectangles functions: Wave Plus only

Blue Rectangles functions: WaveX only

#### WaveX and Waveplus

[DeviceState](#) State { get; }

Represents the current state of the state machine which regularizes the access to the library (see [DeviceState](#))

[DeviceError](#) InitialError { get; }

Includes the codification of a possible error that occurred during the execution of the DaqSystem class construction (see [DeviceError](#))

[List<DeviceType>](#) Type { get; }

Represents the device type list of the WaveX and Waveplus devices connected to USB ports (see [DeviceType](#))

[List<DeviceMode>](#) Mode { get; }

Represents the device mode list of the WaveX and Waveplus devices connected to USB ports (see [DeviceMode](#))

[List<int>](#) DeviceInstalledSensors { get; }

Represents the list of the installed sensors number for each WaveX and Waveplus device connected to USB ports

[List<DeviceDependentFunctionalities>](#) { get; }

It represents the device-dependent functionalities list of the WaveX and Waveplus devices connected to USB ports (see [DeviceDependentFunctionalities](#))

[void](#) ConfigureCapture([ICaptureConfiguration](#) captureConfiguration);

Configures the data acquisition process (see [CaptureConfiguration](#))

[ICaptureConfiguration](#) CaptureConfiguration();

Gets the current configuration of the data acquisition process (see [CaptureConfiguration](#))

[void](#) StartCapturing([DataAvailableEventPeriod](#) dataAvailableEventPeriod);

Determines the start of the data acquisition process  
dataAvailableEventPeriod represents the time interval between two consecutive DataAvailableEvents events (see [DataAvailableEvents](#) and [DataAvailableEventPeriod](#))

[void](#) StopCapturing();

Determines the stop of the data acquisition process

[void](#) GenerateInternalStartTrigger();

Generates via software the internal start trigger activation (see [DataAvailableEventArgs](#))

[void](#) GenerateInternalStopTrigger();

Generates via software the internal stop trigger activation (see [DataAvailableEventArgs](#))

[List<IVersion>](#) FirmwareVersion { get; }

Represents the firmware version list of the WaveX and Waveplus devices connected to the USB ports (see [Version](#))

[List<IVersion>](#) HardwareVersion { get; }

Represents the hardware version list of the WaveX and Waveplus devices connected to the USB ports (see [Version](#))

[IExtVersion](#) SoftwareVersion { get; }

Represents the Waveplus.Daq.NET software library version (see [ExtVersion](#))

**int** InstalledSensors { get; }

Represents the total number of the installed sensors

**int** InstalledFootSwSensors { get; }

Represents the total number of the installed Foot Switch sensors

**SensorModel[]** SensorsModel()

Returns a vector of InstalledSensors elements each of them represents the corresponding sensor model (see [SensorModel](#))

**void** EnableSensor(**int** sensor);

Enables the specified sensor when sensor assumes the value between 1 and InstalledSensors. If sensor assumes the value 0, all the installed sensors are enabled.

**void** DisableSensor(**int** sensor);

Disable (stand-by mode) the selected sensor when sensor assumes the value between 1 and InstalledSensors. If sensor assumes the value 0, all sensors are disabled.

**void** ConfigureSensor([ISensorConfiguration](#) sensorConfiguration, **int** sensor);

Configures the selected sensor when sensor assumes the value between 1 and InstalledSensors. If sensor assumes the value 0, all sensors are configured (see [SensorConfiguration](#))

[ISensorConfiguration](#) SensorConfiguration(**int** sensor);

Returns the current configuration of the selected sensor from the “sensor” parameter. This parameter can have value between 1 and InstalledSensors (see [SensorConfiguration](#))

**event** [EventHandler<DeviceStateChangedEventArgs>](#) StateChanged;

Represents the handler of the StateChanged event. This is generated every time the state machine, that manages the access to the library, changes state (see [DeviceStateChangedEventArgs](#))

**event** [EventHandler<DataAvailableEventArgs>](#) DataAvailable;

Represents the handler of the DataAvailable event. This is generated during the acquisition process, when new samples are available (see Scheme 2). To optimize the data processing and transfer, the samples are buffered and published around every DataAvailableEventPeriod ms (see [StartCapturing\(\)](#) and [DataAvailableEventPeriod](#))

**void** UpdateDisplay();

Update the WaveX display according to the current device configuration

**void** EnableFootSwSensors();

Enables all Foot Switch sensors (2 in total)

**void** DisableFootSwSensors();

Disables (stand-by mode) all Foot Switch sensors (2 in total)

**void** DetectAccelerometerOffset(**int** sensor);

Reads and compensates the offset of x, y, z channels of the selected accelerometers (Waveplus EMG sensor accelerometers), when sensor has value between 1 and InstalledSensors. If sensor has value 0, than all accelerometers are compensated for offset

[SensorCheckReport\[\]](#) CheckElectrodeImpedance(**int** sensor);

Executes the impedance check on the sensor selected, when sensor has value between 1 and InstalledSensors. If sensor has value 0, than all sensors are checked for impedance. It returns a vector whose elements are of the type [SensorCheckReport](#); it includes all check results for all installed sensors (see [SensorCheckReport](#))

**void** CalibrateSensorImuOffset (**int** sensor);

Executes the accelerometer and gyroscope offset calibration of the selected inertial sensor, when sensor has value between 1 and InstalledSensors . If sensor has value 0, than it executes the accelerometer and gyroscope offset calibration of all the installed inertial sensors

**void** CalibrateGyroscopeOffset (**int** sensor);

Executes the gyroscope offset calibration of the selected inertial sensor, when sensor has value between 1 and InstalledSensors . If sensor has value 0, than it executes the gyroscope offset calibration of all the installed inertial sensors

**void** TurnSensorLedOn(**int** sensor);

Activates the blinking of the sensor LED, when sensor has value between 1 and InstalledSensors. If sensor has value 0, than the operation is repeated on all sensors

**void** TurnAllSensorLedsOn();

Activated the blinking of all sensors, including Footswitch sensors

**void** TurnAllSensorLedsOff();

Disable the blinking of all sensors, including Footswitch sensors


## 3.2 DeviceState

Defines the ensemble of all states that the state machine can assume:

**public enum** DeviceState

{	
NotConnected,	No USB device is connected to the bus USB
Initializing,	WaveX /Waveplus is initializing
CommunicationError,	No communication with WaveX/Waveplus device available
InitializingError,	WaveX/Waveplus device was not correctly initialized
Idle,	WaveX/Waveplus device is ready for use
Capturing,	WaveX/Waveplus device is acquiring data
ReadingSensorMemory,	WaveX/Waveplus device is reading data from sensor
memory	
}	

Typical state transition and events generation during data capturing process



### 3.3 CaptureConfiguration

The class [CaptureConfiguration](#) implements the interface [ICaptureConfiguration](#) that includes the following properties:

[bool](#) ExternalTriggerEnabled { [get](#); [set](#); }  
Enables/disables the external trigger

[int](#) ExternalTriggerActiveLevel { [get](#); [set](#); }  
Defines the level at which trigger will start/stop the acquisition

[bool](#) ExternalTriggerStopDisable { [get](#); [set](#); }  
Disables/Enables the external stop trigger

[ImuAcqXType](#) IMU\_AcqXType { [get](#); [set](#); }  
Defines WaveX IMU sensors acquisition type (see [ImuAcqXType](#))

[EmgAcqXType](#) EMG\_AcqXType { [get](#); [set](#); }  
Defines WaveX EMG sensors acquisition type (see [EmgAcqXType](#))

[EmgImuAcqXType](#) EMG\_IMU\_AcqXType { [get](#); [set](#); }  
Defines WaveX EMG + IMU sensors acquisition type (see [EmgImuAcqXType](#))

[IFootSwTransducerEnabled](#) FootSwATransducerEnabled { [get](#); [set](#); }  
Enables/Disables the Footswitch A transducers (see [FootSwTransducerEnabled](#))

[IFootSwTransducerEnabled](#) FootSwBTransducerEnabled { [get](#); [set](#); }  
Enables/Disables the Footswitch B transducers (see [FootSwTransducerEnabled](#))

[IFootSwTransducerThreshold](#) FootSwATransducerThreshold { [get](#); [set](#); }  
Defines the activation levels for the transducers of Footswitch A (see [FootSwTransducerThreshold](#))

[IFootSwTransducerThreshold](#) FootSwBTransducerThreshold { [get](#); [set](#); }  
Defines the activation levels for the transducers of Footswitch B (see [FootSwTransducerThreshold](#))

[FootSwProtocol](#) FootSwProtocol { [get](#); [set](#); }  
Defines the protocol of processing of the data coming from the Footswitch sensors  
(see [FootSwProtocol](#))

[ImuAcqPType](#) IMU\_AcqPType { [get](#); [set](#); }  
Defines Waveplus IMU sensors acquisition type (see [ImuAcqPType](#))

[EmgAcqPType](#) EMG\_AcqPType { [get](#); [set](#); }  
Defines Waveplus EMG sensors acquisition type (see [EmgAcqPType](#))

Defines the WaveX inertial sensors acquisition type

[enum](#) [ImuAcqXType](#)

```
{
    RawAccGyroMagData_400Hz    raw data acquisition (accelerometer, gyroscope, magnetometer) at
                                400 Hz
    RawAccGyroData_500Hz       raw data acquisition (accelerometer, gyroscope) at 500 Hz
    Fused6xData_400Hz          6 axis fused data (quaternions) acquisition at 400 Hz
    Fused9xData_400Hz          9 axis fused data (quaternions) acquisition at 400 Hz
    Mixed6xData_200Hz          6 axis fused data (quaternions) acquisition at 200 Hz
                                and raw data acquisition (accelerometer, gyroscope at 200 Hz,
                                magnetometer at 100 Hz)
    Mixed9xData_200Hz          6 axis fused data (quaternions) acquisition at 200 Hz
                                and raw data acquisition (accelerometer, gyroscope at 200 Hz,
                                magnetometer at 100 Hz)
}
```

**EmgAcqXType**

Defines the WaveX EMG and accelerometer sensors acquisition type

[enum EmgAcqPType](#)

```
{
    Emg_2kHz                EMG acquisition at 2 kHz
}
```

**EmgImuAcqXType**

Defines the WaveX EMG+IMU sensors acquisition type

[enum EmgImuAcqXType](#)

```
{
    Emg_2kHz                EMG acquisition at 2 kHz
    Emg_2kHz_RawAccGyroMagData_100Hz EMG acquisition at 2 kHz and raw data acquisition
                                   (accelerometer, gyroscope and magnetometer at 100 Hz)
    Emg_2kHz_RawAccGyroData_200Hz    EMG acquisition at 2 kHz and raw data acquisition
                                   (accelerometer and gyroscope at 200 Hz)
    Emg_2kHz_Fused6xData_100Hz       EMG acquisition at 2 kHz and 6 axis fused data
                                   (quaternions) acquisition at 100 Hz
    RawAccGyroMagData_400Hz           raw data acquisition (accelerometer, gyroscope,
                                   magnetometer) at 400 Hz
    RawAccGyroData_500Hz              raw data acquisition (accelerometer, gyroscope) at 500 Hz
    Fused6xData_400Hz                6 axis fused data (quaternions) acquisition at 400 Hz
    Fused9xData_400Hz                9 axis fused data (quaternions) acquisition at 400 Hz
    Mixed6xData_200Hz                6 axis fused data (quaternions) acquisition at 200 Hz
                                   and raw data acquisition (accelerometer, gyroscope at 200
                                   Hz, magnetometer at 100 Hz)
    Mixed9xData_200Hz                6 axis fused data (quaternions) acquisition at 200 Hz
                                   and raw data acquisition (accelerometer, gyroscope at 200
                                   Hz, magnetometer at 100 Hz)
}
```

8

**ImuAcqPType**

Defines the Waveplus inertial sensors acquisition type

[enum ImuAcqXType](#)

```
{
    RawAccGyroMagData_284Hz    raw data acquisition at 284 Hz
    Fused9xData_142Hz          9 axis fused data (quaternions) acquisition at 142 Hz
    Fused6xData_284Hz          6 axis fused data (quaternions) acquisition at 284 Hz
    Fused9xData_71Hz           9 axis fused data (quaternions) acquisition at 71 Hz
    Fused6xData_142Hz          6 axis fused data (quaternions) acquisition at 142 Hz
    Mixed6xData_142Hz          6 axis fused data (quaternions) acquisition at 142 Hz
                                   and raw data acquisition (accelerometer and gyroscope at 142 Hz,
                                   magnetometer at 47 Hz)
}
```

**EmgAcqPType**

Defines the Waveplus EMG and accelerometer sensors acquisition type

[enum EmgAcqPType](#)

```
{
    Emg_2kHz_Acc_142Hz        EMG acquisition at 2 kHz, accelerometer acquisition at 142 Hz
}
```



### FootSwTransducerEnabled

The class [FootSwTransducersEnabled](#) implements the interface [IFootSwTransducersEnabled](#) that includes the following properties:

```
bool T_A { get; set; }
    Enables/Disables the transducer A

bool T_1 { get; set; }
    Enables/Disables the transducer 1

bool T_5 { get; set; }
    Enables/Disables the transducer 5

bool T_T { get; set; }
    Enables/Disables the transducer T
```

### FootSwTransducerThreshold

The class [FootSwTransducersThreshold](#) implements the interface [IFootSwTransducersThreshold](#) that includes the following properties:

```
double T_A { get; set; }
    Defines the threshold used for transducer A

double T_1 { get; set; }
    Defines the threshold used for transducer 1

double T_5 { get; set; }
    Defines the threshold used for transducer 5

double T_T { get; set; }
    Defines the threshold used for transducer T
```

### FootSwProtocol

Defines the protocol used to process the Footswitch samples

```
public enum FootSwProtocol
{
    FullFoot,      FullFoot protocol
    HalfFoot,      HalfFoot protocol
    QuarterFoot    QuarterFoot protocol
}
```

9

## 3.9 Version

The class [Version](#) implements the interface [IVersion](#) that includes the following properties:

```
int Major { get; }
    Represents the Major part of the version number

int Middle { get; }
    Represents the Middle part of the version number

int Minor { get; }
    Represents the Minor part of the version number
```

## 3.10 ExtVersion

The class [ExtVersion](#) implements the interface [IExtVersion](#) that includes the following properties:

```
int Major { get; }
    Represents the Major part of the version number
```

```
int Minor { get; }
```

Represents the Minor part of the version number

```
int Build { get; }
```

Represents the Build part of the version number

```
int Revision { get; }
```

Represents the Revision part of the version number

### 3.11 SensorConfiguration

The class [SensorConfiguration](#) implements the interface [ISensorConfiguration](#) that includes the following properties:

```
SensorModel SensorModel { get; set; }
```

Represents the sensor mode (see [SensorModel](#))

```
SensorMode SensorMode { get; set; }
```

Represents the sensor mode (see [SensorMode](#))

```
AccelerometerFullScale AccelerometerFullScale { get; set; }
```

Represents the full scale value for the axis x, y, z of the accelerometer (see [AccelerometerFullScale](#))

```
GyroscopeFullScale GyroscopeFullScale { get; set; }
```

Represents the full scale value for the axis x, y, z of the gyroscope (see [GyroscopeFullScale](#))

### 3.12 SensorModel

Represents the sensor model:

```
enum SensorModel
```

```
{
```

Undefined	undefined
Mini_EmgImu	Mini EMG + IMU
Mini_Emg	Mini EMG
Pico_EmgImu	Pico EMG + IMU
Pico_Emg	Pico EMG
Imu	Imu

MiniPicoImuWavePlus	Waveplus sensor model
---------------------	-----------------------

```
}
```

10

### 3.12 SensorMode

Represents the sensor mode:

```
enum SensorMode
```

```
{
```

EMG_SENSOR	EMG
INERTIAL_SENSOR	inertial
ANALOG_GP_SENSOR	general purpose

EMG_INERTIAL_SENSOR	EMG + inertial
---------------------	----------------

FSW_SENSOR	foot-switch
EMG_ACC_SENSOR	EMG + accelerometer

```
}
```

### 3.13 AccelerometerFullScale

Represents full scale value for the axis x, y, z of the accelerometer

```
public enum AccelerometerFullScale
{
    g_2,           Full scale = 2g
    g_4,           Full scale = 4g
    g_8,           Full scale = 8g
    g_16          Full scale = 16g
}
```

g = gravitational field

### 3.14 GyroscopeFullScale

Represents full scale value for the axis x, y, z of the gyroscope

```
public enum GyroscopeFullScale
{
    dps_250,       Full scale = 250 D/s
    dps_500,       Full scale = 500 D/s
    dps_1000,      Full scale = 1000 D/s
    dps_2000       Full scale = 2000 D/s
}
```

D/s = degree per second

### 3.15 SensorCheckReport

Defines the result of the impedance check of one sensor

```
public enum SensorCheckReport
{
    NotPassed,
    Passed,
    NotExecuted           (the sensor was not enabled or the test was not requested for that sensor)
}
```

11

### 3.16 DeviceStateChangedEventArgs

The class [DeviceStateChangedEventArgs](#) implements the following property:

```
DeviceState State { get; }
    Represents the current device state
```

### 3.16 DataAvailableEventArgs

The class [DataAvailableEventArgs](#) implements the following properties:

```
public int ScanNumber { get; set; }
    Represents the number of samples available for each channel

public float[,] EmgSamples;
    Represents a vector with two dimensions including ScanNumber samples for each EMG channel
    The samples are expressed in [uV] unit
    The first index identifies the sensor (it can be between 0 and InstalledSensors -1)
    The second index identifies the sample (it can be between 0 and ScanNumber -1)

public float[,] ImuSamples;
    Represents a vector with three dimensions including ScanNumber samples for each IMU quaternion
    components
    The first index identifies the sensor (it can be between 0 and InstalledSensors -1)
    The second can assume values of 0, 1, 2, 3 that identifies respectively the components
    w, x, y and z of the quaternion
    The third index identifies the sample (it can be between 0 and ScanNumber -1)
    ImuSamples samples are available only during IMU fused data acquisition or mixed acquisition

public float[,] AccelerometerSamples;
```

Represents a vector with three dimensions including ScanNumber samples for each Accelerometer channel  
 The samples are expressed in [g] unit  
 The first index identifies the sensor (it can be between 0 and InstalledSensors -1)  
 The second can assume values of 0, 1, 2 that identifies respectively the values x, y and z of the acceleration vector  
 The third index identifies the sample (it can be between 0 and ScanNumber -1)

**public float**[,] GyroscopeSamples;

Represents a vector with three dimensions including ScanNumber samples for each Gyroscope channel  
 The samples are expressed in [D/s] unit  
 The first index identifies the sensor (it can be between 0 and InstalledSensors -1)  
 The second can assume values of 0, 1, 2 that identifies respectively the values x, y and z of the angular velocity vector  
 The third index identifies the sample (it can be between 0 and ScanNumber -1)

**public float**[,] MagnetometerSamples;

Represents a vector with three dimensions including ScanNumber samples for each Magnetometer channel  
 The samples are expressed in [uT] unit  
 The first index identifies the sensor (it can be between 0 and InstalledSensors -1)  
 The second can assume values of 0, 1, 2 that identifies respectively the values x, y and z of the magnetic vector  
 The third index identifies the sample (it can be between 0 and ScanNumber -1)

**public short**[,] FootSwSamples;

Represents a vector with two dimensions including the ScanNumber for each Footswitch channel  
 The samples are expressed in [V] unit  
 The first number identifies the sensor (it can be between 0 and InstalledFSWSensors -1)  
 The second number identifies the sample (it can be between 0 and ScanNumber -1)

**public float**[] SyncSamples;

Reserved

12

**public short**[] SensorStates;

Represents a vector including the state of each sensor .  
 The index value identifies the sensor (it can be between 0 and InstalledSensors -1)  
 (See Sensors State)

**public short**[] FootSwSensorStates;

Represents a vector including the state of each Foot-switch  
 The index value identifies the sensor (it can be between 0 and InstalledFSWSensors -1)  
 (See Sensors State)

**public bool** StartTriggerDetected;

Indicates if a start trigger signal was detected

**public bool** StopTriggerDetected;

Indicates if a stop trigger signal was detected

**public int** StartTriggerScan;

Represents the position (sample index value), in the data vectors, of the sample in correspondence of which the start trigger was detected. If no start trigger was detected, the value is 0

**public int** StopTriggerScan;

Represents the position (sample index value), in the data vectors, of the sample in correspondence of which the stop trigger was detected. If no stop trigger was detected, the value is 0

### 3.17 DataAvailableEventPeriod

Represents the time interval between two consecutive DataAvailableEvents events

```

public enum DataAvailableEventPeriod
{
    ms_100 ,      Interval = 100 ms
    ms_50 ,       Interval = 50 ms
    ms_25 ,       Interval = 25 ms
    ms_10        Interval = 10 ms
}

```

**Note:** interval values less than 100 ms need a very short execution-time of DataAvailableEvent event handler

### 3.18 Sensors state

Sensor state is represented by a 16 bit formatted according to the following table:

D15	D14	D13	D12	D11	D10	D9	D8	D7	D6	D5	D4	D3	D2	D1	D0
0	0	0	0	0	0	0	0	0	0	0	BLAF	BL3	BL2	BL1	BL0

where BL3 – BL0 represent the sensors battery charge level, BLAF represents the battery level available flag:

BLA	BL3	BL2	BL1	BL0	Battery charge level
1	0	0	0	0	0%
1	0	0	0	1	
1	0	0	1	0	
1	0	0	1	1	
1	0	1	0	0	
1	0	1	0	1	
1	0	1	1	0	
1	0	1	1	1	
1	1	0	0	0	
1	1	0	0	1	
1	1	0	1	0	
1	1	0	1	1	
1	1	1	0	0	
1	1	1	0	1	
1	1	1	1	0	
1	1	1	1	1	100%
0	X	X	X	X	Not Available

**Note:** for Waveplus system the sensors battery charge level information is available only during EMG/Accelerometer, Foot-switch and IMU raw-data acquisitions

### 3.19 DeviceError

Defines the possible errors:

```

enum DeviceError
{
    Success = 0,
    UnspecifiedError= 1,
    DeviceNotConnected = 2,
    SendingCommand = 3,
    ReceivingCommandReply = 4,
    DeviceErrorExecutingCommand = 5,
    ConfiguringCapture = 6,
    ConfiguringAudio = 7,
    WrongCaptureConfigurationImuAcqType = 8,
    WrongCaptureConfigurationSamplingRate = 9,
    WrongCaptureConfigurationSamplingRateFromDevice = 10,
    WrongCaptureConfigurationDataAvailableEventPeriod = 11,
    WrongAudioConfigurationSpeakerLevel = 12,
    WrongAudioConfigurationSource = 13,
}

```

WrongAudioConfigurationSourceFromDevice = 14,  
WrongSensorConfigurationAccelerometerFullScale = 15,  
WrongSensorConfigurationGyroscopeFullScale = 16,  
WrongSensorConfigurationSensorType = 17,  
WrongFootSwProtocol = 18,  
WrongDeviceTypeFromDevice = 19,  
WrongDeviceMode = 20,  
WrongSensorNumber = 21,  
WrongFootSwSensorNumber = 22,  
FootSwSensorNotInstalled = 23,  
ReadingBackCaptureConfigurationSamplingRate = 24,  
ReadingBackAudioConfigurationAudioSource = 25,  
ReadingBackAudioConfigurationSpeakerVolumeLevel = 26,  
ReadingBackAudioConfigurationSpeakerEnabled = 27,  
ReadingBackCommunicationTestData = 28,  
ReadingBackSensorCommandBuffer = 29,  
ParsingFpgaConfigurationFile = 30,  
SendingFpgaConfigurationImage = 31,  
ReceivingUpdateFpgaConfigurationCommandReply = 32,  
DeviceErrorExecutingUpdateFpgaConfigurationCommand = 33,  
FpgaConfiguredFlagNotAsserted = 34,  
ReadingFpgaConfiguredFlag = 35,  
DataTransferThreadStartingTimeout = 36,  
ActionNotAllowedInTheCurrentDeviceState = 37,  
FeatureNotSupportedByTheDevice = 38,  
ActionNotAllowedInTheCurrentDataTransferState = 39,  
WrongTrialNumber = 40,  
WrongDeviceState = 41,  
WrongDeviceAction = 42,  
TimeoutExecutingSensorCommand = 43,  
CommandNotExecutedByAllTheSensors = 44,  
TimeoutReadingSensorMemoryStatus = 45,  
SensorMemoryHeaderCorrupted = 46,  
SensorMemoryWrongTrialAddressFormat = 47,  
SensorMemoryTrialAddressOutOfRange = 48,  
WrongDaqTimeOutValue = 49,  
BadSensorCommunication = 50,  
SyncBuffer1Overrun = 51,  
SyncBuffer2Overrun = 52,  
ImuCalibrationNotAvailable = 53,  
SaveMPStatusNotAvailable = 54,  
FunctionNotAvailable = 55,  
WrongDeviceRFChannelFromDevice = 56,  
WrongImuAcqType = 57,  
SensorDoesNotProvideSelectiveMemoryReading = 58,  
ConflictingDeviceModeFromEEPROMandFromMicro = 59,  
WrongEmgHwGain = 60,  
WrongCaptureConfigurationImuAcqPType = 61,  
WrongCaptureConfigurationImuAcqXType = 62,  
WrongCaptureConfigurationEmgAcqPType = 63,  
WrongCaptureConfigurationEmgAcqXType = 64,  
WrongCaptureConfigurationEmgImuAcqXType = 65,  
WrongCaptureConfigurationEmgAcqBType = 66,  
WrongCaptureConfigurationImuAcqBType = 67,  
UnavailableBSystemEmgAcqType = 68,  
UnavailablePSystemEmgAcqType = 69,  
UnavailablePSystemImuAcqType = 70,  
  
SendingAnalogOutputsConfiguration = 71,  
ReceivingAnalogOutputsConfigurationReply = 72,  
DeviceErrorSendingAnalogOutputsConfiguration = 73,  
ReadingBackAnalogOutputsConfiguration = 74,  
  
SendingSensorMemoryStatusRequest = 75,  
ReceivingSensorMemoryStatusReply = 76,

```

ExecutingSensorMemoryStatusRequest = 77,
DataPacketLost = 78,

WiFiSensorMemoryReadingError1 = 79,
WiFiSensorMemoryReadingError2 = 80,

SendingSensorIdRequest = 81,
ReceivingSensorIdReply = 82,
ExecutingSensorIdRequest = 83,

WrongSensorConfigurationSensorId = 84,

WrongRFDeviceChannel = 85,
WrongRFSensorChannel = 86,

WiFiSensorMemoryReadingError3 = 87,
WiFiSensorMemoryReadingError4 = 88,
WiFiSensorMemoryReadingError5 = 89,
WiFiSensorMemoryReadingError6 = 90,
WiFiSensorMemoryReadingError7 = 91,
WiFiSensorMemoryReadingError8 = 92,

WrongDeviceCurrentTrialIdFromDevice = 93,
WrongCommandCodeFromSensor = 94,
WrongRFChannel = 95,
WrongRFChannelName = 96,
WrongSensorId = 97,
WrongSensorIdCode = 98,
WrongSensorIdName = 99,
}

```

### 3.20 DaqDeviceExceptionType

15

Defines the possible exception types generated by the library:

```

enum DaqDeviceExceptionType
{
    deviceNotConnected = 0,
    unableToStartCaptureDataTransfer = 1,
    unableToStartImpedanceDataTransfer = 2,
    unableToStartCapturing = 3,
    unableToStopCapturing = 4,
    unableToReadSensorMemoryStatus = 5,
    unableToGetCaptureConfiguration = 6,
    unableToSetCaptureConfiguration = 7,
    unableToGetInstalledSensors = 8,
    unableToGetDeviceType = 9,
    unableToConfigureSensor = 10,
    unableToGetSensorConfiguration = 11,
    unableToTurnInternalTrigger_OFF = 12,
    unableToTurnInternalTrigger_ON = 13,
    unableToEnableSensor = 14,
    unableToDisableSensor = 15,
    unableToEnableFootSwSensor = 16,
    unableToDisableFootSwSensor = 17,
    unableToDetectAccelerometerOffset = 18,
    unableToCheckElectrodeImpedance = 19,
    unableToGetElectrodeImpedanceReport = 20,
    unableToTurnSensorLedOn = 21,
    unableToTurnFootSwSensorLedOn = 22,
    unableToTurnAllSensorLedsOn = 23,
    unableToTurnAllSensorLedsOff = 24,
    unableToStartSensorMemoryRecording = 25,
    unableToStopSensorMemoryRecording = 26,
}

```

```

unableToClearSensorMemory = 27,
unableToFormatSensorMemory = 28,
unableToStartSensorMemoryReading = 29,
unableToStartSensorSelectiveMemoryReading = 30,
unableToStopSensorMemoryReading = 31,
unableToEnableSensorMemoryMode = 32,
unableToDisableSensorMemoryMode = 33,
unableToCalibrateSensorImu = 34,
unableToCalibrateSensorGyroscope = 35,
unableToSaveMPStatus = 36,
unableToGetFirmwareVersion = 37,
unableToGetHardwareVersion = 38,
unableToSetAudioConfiguration = 39,
unableToGetAudioConfiguration = 40,
unableToConvertParameter = 41,
unableToUpdateFirmware = 42,
unableToGetFPGAConfigFlag = 43,
unableToUpdateDeviceBoardEEPROMInfo = 44,
unableToGetDeviceBoardEEPROMInfo = 45,
unableToSynchronizeData = 46,
unableToChangeDeviceRFChannel = 47,
unableToChangeSensorRFChannel = 48,
unableToGetDeviceRFChannel = 49,
unableToSetFirstImuCalibrationStep = 50,
unableToSetNextImuCalibrationStep = 51,
unableToSetIMUAcqType = 52,
unableToGetDeviceDependentFunctAvailability = 53,
unableToReadSensorInfo = 54,
unableToEnableRemoteControl = 55,
unableToDisableRemoteControl = 56,
unableToCalibrateSensorAccelerometer = 57,
unableToSetEmgImuAcqType = 58,
unableToGetDeviceSensorsId = 59,
unableToSetEmgAcqType = 60,
unableToSetImuAcqType = 61,
unableToSetAnalogOutputsConfiguration = 62,
unableToGetAnalogOutputsConfiguration = 63,
unableToEnableDeviceWiFi = 64,
unableToDisableDeviceWiFi = 65,
unableToEnableSensorWiFi = 66,
unableToDisableSensorWiFi = 67,
unableToSetSensorMemoryTrialId = 68,
unableToStartSensorMemoryDataTransfer = 69,
unableToFindAndInitializeBDevice = 70,
unableToDisconnectBDevice = 71,
unableToSetDeviceDmaBufferSize = 72,
unableToDisableWiFi = 73,
unableToGetSensorMemoryStatus = 74,
unableToReadSensorId = 75,
unableToSetSensorSleepMode = 76,
unableToSaveConfigurationOnXDevice = 77,
unableToGetDeviceCurrentTrialId = 78,
unableToGetSensorInfo = 79,
unableToCalibrateSensorEmg = 80,
unableToGetFixedMemoryId = 81,
}

```

16

#### 4. State Machine: transitions and availability of Properties and Methods

The access to the library functionalities is regulated by a state machine aimed at improving the software reliability. In Scheme 3 the allowed transitions and the available methods/properties are highlighted for every state.



Scheme 3

